


TOPICAL SKIN ADHESIVES

S2S Global has partnered with Chemence Medical, an established global medical device market leader, to offer an expansive product line of PremierPro™ Topical Skin Adhesives.

PremierPro™ EXOFIN® Topical Skin Adhesive (TSA) is a sterile tissue adhesive utilizing a fast-drying monomeric formulation. The high viscosity product reduces the risk of unintended placement of the adhesive, and the colorant contained in the formulation allows the user to easily see where the adhesive has been applied. Sold in 1.0 ml and 0.5 ml Micro versions.

PremierPro™ EXOFIN FUSION® is a complete sterilized Skin Closure System for longer incision wounds.



OFFERINGS

 PremierPro™
exofin®

NEW!

 PremierPro™
exofin® PRECISION
PEN PremierPro™
exofin®
MICRO PremierPro™
exofin fusion®
SKIN CLOSURE SYSTEM

EXOFIN® and EXOFIN FUSION® are registered trademarks of Chemence Medical, Inc.
Please visit chemencemedical.com for source data regarding all product performance claims.



TOPICAL SKIN ADHESIVE (TSA)

FAST

Fastest polymerizing 2-octyl topical adhesive on the market today. Average set times are less than 30 seconds.

STRONG

2-Octyl formulation is proven to be the strongest and most flexible adhesive available. Creates a microbial barrier against clinically relevant pathogens. Added flexibility allows the adhesive to remain intact for 7-10 days.

INNOVATIVE

Soft and flexible applicator glides across uneven surfaces, creating a solid seal on skin. Transparent tip allows user to see the adhesive prior to expression.

EASY TO USE

Simply screw on the pre-assembled applicator tip and begin to express adhesive. Apply as much or as little adhesive as needed, and get the same great results every time!

VALUE

1 ml tube fill offers 25% more adhesive than the largest adhesive available on the market*. Decreases unit volume, thus adding cost savings per procedure.

*data available upon request





DESCRIPTION	SIZE	REORDER NUMBER	PACKAGING
High Viscosity Tissue Adhesive	1.0 ml	3470	10 tubes/bx

HIGH VISCOSITY TOPICAL SKIN ADHESIVE PEN



- + Industry leading **14 day microbial barrier**
- + **Highest viscosity** in industry
- + **Fastest drying adhesive** in industry
- + **Most usable adhesive** in industry (152% more)
- + Intuitive, point-and-shoot design with multiple application techniques
- + See-through activation chamber
- + Ergonomic design with soft non-slip grip for precise control
- + Pressure controlled application for metered dose

THREE APPLICATION OPTIONS


Broad Coverage



Narrow LineP



Precise Line



DESCRIPTION	SIZE	REORDER NUMBER	PACKAGING
High Viscosity Topical Skin Adhesive Pen	1.0 ml	3475	12 pens/bx



TOPICAL SKIN ADHESIVE (TSA) - MICRO

- + Ideal for **small trauma induced injuries** and **single port incisions**.
- + **Innovative silicone tip** perfectly suited for precision application.
- + **Clear applicator** allows user to see the adhesive prior to expression from tip.
- + **Same high viscosity formulation** as full size PremierPro™ EXOFIN® Topical Skin Adhesive (TSA).



DESCRIPTION	SIZE	REORDER NUMBER	PACKAGING
High Viscosity Tissue Adhesive (Micro)	0.5 ml	3471	10 tubes/bx

SKIN CLOSURE SYSTEMS



- + EXOFIN FUSION® is an innovative two-part system intended for use in three different sizes for **medium to long incisions up to 58 cm in length***.
- + Products contain **polyester mesh tape and 2-octyl cyanoacrylate monomer**.
- + Typical adhesive **cure time is less than one minute** compared to average closure using sutures of 7 to 20 minutes depending upon procedure.
- + Provides **microbial barrier** preventing bacteria from entering incision for up to 21 days*.
- + **Distributes tension along the entire incision** allowing uniformity along the line, thus preventing skin gaps from forming.
- + Curved ends **reduce the area for potential skin damage** due to stress and strain by 67%, dissipating friction and force across the entirety of the strip.
- + No initiators in the mesh itself which **reduces patient exposure to initiators** by 61%. Instead, the initiator is delivered with the polymerizable adhesive.

*Data available upon request

Each Skin Closure System Contains:

22cm (Item 3472)

- (1) 4cm x 22cm Mesh Patch
- (2) 1.75g Tubes of Adhesive

30cm (Item 3474)

- (1) 4cm x 30cm Mesh Patch
- (2) 1.75g Tubes of Adhesive

60cm (Item 3473)

- (2) 4cm x 30cm Mesh Patches
- (3) 1.75g Tubes of Adhesive



DESCRIPTION	SIZE	REORDER NUMBER	PACKAGING
Skin Closure System (2 per box): (1) 4cm x 22cm Mesh Patch (2) 1.75g Tubes of Adhesive	22 cm	3472	2 kits/bx
Skin Closure System (2 per box): (1) 4cm x 30cm Mesh Patch (2) 1.75g Tubes of Adhesive	30 cm	34742	2 kits/bx
Skin Closure System (2 per box): (2) 4cm x 30cm Mesh Patches (3) 1.75g Tubes of Adhesive	60 cm	3473	2 kits/bx