



Heel Incision Safety Lancets



New integrated safety tab with locking trigger for one hand operation




Sculpted design with anti-slip texture for easier grip



Tapered lancing position with embossed arrow to improve visibility and accuracy


DESCRIPTION	COLOR	REFERENCE NUMBER	PACKAGING
Heel Incision Safety Lancets Depth: 0.65mm/Length: 1.4mm	Green	7511	1/pk; 50/bx; 4 bxs/cs
Heel Incision Safety Lancets Depth: 0.85mm/Length: 1.75mm	Orange	7512	1/pk; 50/bx; 4 bxs/cs
Heel Incision Safety Lancets Depth: 1mm/Length: 2.5mm	Yellow	7513	1/pk; 50/bx; 4 bxs/cs

Pressure Activated Safety Lancets




DESCRIPTION	COLOR	REFERENCE NUMBER	PACKAGING
Pressure Activated Safety Lancets 21G Needle x 2.0mm Depth	Pink	7551	100/bx; 20 bxs/cs
Pressure Activated Safety Lancets 25G Needle x 1.8mm Depth	Blue	7552	100/bx; 20 bxs/cs
Pressure Activated Safety Lancets 30G Blade x 1.5mm Depth	Light Green	7553	100/bx; 20 bxs/cs

Instructions for Use




Button Activated Safety Lancets



DESCRIPTION	COLOR	REFERENCE NUMBER	PACKAGING
Button Activated Safety Lancets 17G Needle x 2.0mm Depth	Yellow	7571	100/bx; 20 bxs/cs
Button Activated Safety Lancets 23G Needle x 1.8mm Depth	Blue	7572	100/bx; 20 bxs/cs
Button Activated Safety Lancets 28G Blade x 1.5mm Depth	Purple	7573	100/bx; 20 bxs/cs

Instructions for Use



Twist off the protective cap and pull it straight out.



Place Lancet on the desired puncture site and press the top push button to activate the device.



Gently apply intermittent pressure near the puncture site to obtain the required blood volume.