

BlueLine™


Premium Repositioning Underpad



EXPECT MORE in Moisture Management

VERSATILE


- All-in-one solution for moisture management and patients' turns, boosts, and repositioning
- Air-permeable backsheet is for use on low air loss mattresses and specialty support surfaces
- Exceeds National Pressure Injury Advisory Panel (NPIAP) recommendation for breathability
- Fight skin breakdown with powerful absorbing technology
- Helps prevent pressure injuries, incontinence associated dermatitis (IAD), and moisture associated skin damage (MASD)



Advanced superabsorbent core and diamond embossment design both work together to wick fluid away from patient's skin

Soft top sheet provides comfort to patient

Air permeable allowing air to circulate



High-strength breathable backsheet to help assist with patient repositioning

*Image of SKU 676111

* Product images are provided for illustrative purposes only. Due to ongoing branding and packaging updates, product images may depict legacy PremierPro® branding or other packaging variations and may not represent the actual product labeling, artwork, or branding currently in distribution. Product specifications remain as stated unless otherwise noted.


EFFECTIVE

- Offers high absorption and retention capacity while maintaining an optimal microclimate¹
- Easy to reposition patients on a premium underpad with a strength claim which eliminates the need for a drawsheet and reusable underpads
- Superabsorbent polymer (SAP) locks in fluids from incontinent episodes and unpleasant odors in facility
- Custom packaging configuration of 5/BG enables up to 83% waste reduction which can generate savings²

PATIENT AND CLINICIAN FRIENDLY

- Help improve patients' comfort and dignity
- Top sheet feels soft and smooth against skin compared to layers of standard underpad
- Help improve clinicians' efficiency and confidence in facility
- Complements updated infection control protocols in facility³

HOW IT WORKS



Nonwoven Outer Layer

Provides soft cloth-like feels to the underpad and increases strength for repositioning




Superabsorbent Polymer Beads

Absorbs and traps fluid securely which offer a dry feeling for the patient



Wood Pulp Fluff

Offers quick and uniform distribution of fluid into the absorbent core for optimal moisture management



Air-permeable Poly Layer

Prevents fluid from moving through the underpad while allowing air in order to be breathable

REFERENCE NUMBER	DESCRIPTION	SIZE AND RECOMMENDED USE	ABSORBENT CORE	COLOR	PACKAGING
676110	Premium Quilted Ultra Underpad – rated 375 lbs. to reposition patient	23 IN. X 36 IN. <i>Overnight incontinence, moisture protection, postpartum</i>	AIRLAID FLUFF + SAP	BLUE	5/BG, 14 BG/CS
676111	Premium Quilted Ultra Underpad – rated 375 lbs. to reposition patient	30 IN. X 36 IN. <i>Overnight incontinence, moisture protection, postpartum</i>	AIRLAID FLUFF + SAP	BLUE	5/BG, 8 BG/CS
3771	Premium Quilted Ultra Underpad – rated 500 lbs. to reposition patient	40 IN. X 59 IN. <i>Bariatric and EMS/ambulatory patients, used with patient repositioning system (slider sheet)</i>	AIRLAID FLUFF + SAP	WHITE	5/BG, 5 BG/CS

[1]. Data on file. Available upon request

[2]. Data from a collaborative project with an incontinence focus group consisting of Premier members and S2S Global

[3]. Centers for Medicare and Medicaid Services. State Operation Manual # 100-07, Appendix PP, Interpretive Guidelines for Long-Term Care Facilities, 483.80 Infection Control