


BLOOD TRANSFER PRODUCTS



- Easily and securely connects collection tubes to a syringe or sampling port
- Blue Luer Lock (7580) - male connection for use with sampling ports
- Pink Blood Transfer (7581) - female connection for use with syringes
- Not made with natural rubber latex

DESCRIPTION	REFERENCE NUMBER	PACKAGING
Luer Lock Access Device	7580	200/bx; 4 bxs/cs
Blood Transfer Device	7581	200/bx; 4 bxs/cs

SAFETY LANCETS



- High, medium, or low flow lancets available for all types of capillary blood sampling
- Preloaded design with tri-bevel needle or bi-level blade for high speed penetration
- Needle retracts automatically after use to reduce risk of cross contamination
- Color coding helps health care professionals choose the correct lancet and keep lancets organized
- 17G blade option available in button activated lancets
- All lancets are sterile and disposable for single use
- Not made with natural rubber latex

DESCRIPTION	COLOR	REFERENCE NUMBER	PACKAGING
Heel Incision Safety Lancets Depth: 0.65mm/Length: 1.4mm/Micro-Preemie	Green	7511	1/pk; 50/bx; 4 bxs/cs
Heel Incision Safety Lancets Depth: 0.85mm/Length: 1.75mm/Preemie	Orange	7512	1/pk; 50/bx; 4 bxs/cs
Heel Incision Safety Lancets Depth: 1.0mm/Length: 2.5mm/Newborn	Yellow	7513	1/pk; 50/bx; 4 bxs/cs
Pressure Activated Safety Lancets 21G Needle x 2.0mm Depth	Pink	7551	100/bx; 20 bxs/cs
Pressure Activated Safety Lancets 25G Needle x 1.8mm Depth	Blue	7552	100/bx; 20 bxs/cs
Pressure Activated Safety Lancets 30G Needle x 1.5mm Depth	Light Green	7553	100/bx; 20 bxs/cs
Button Activated Safety Lancets 17G Blade x 2.0mm Depth	Yellow	7571	100/bx; 20 bxs/cs
Button Activated Safety Lancets 23G Needle x 1.8mm Depth	Blue	7572	100/bx; 20 bxs/cs
Button Activated Safety Lancets 28G Needle x 1.5mm Depth	Purple	7573	100/bx; 20 bxs/cs